


## SCC-LIS\_TIP OF THE WEEK

MODULE: SoftLab/SoftMic	SOFTWARE VERSION: 4.0.8						
TOPIC: Restoring the Diffpad	ISSUE DATE: 20Jun2018						
DISTRIBUTION: Technologists	PAGE: 1 of 1						
ISSUED BY: M-L Dufresne, Territorial LIS Administrator							

## **TOPIC/QUESTION:**

My Diffpad is gone and all the Diffpad menu options are greyed out?!?!?!?!

How do I get my Diffpad back!?!?



## ANSWER/TIP:

Click the "On\Off Diffpad" Button



## EXAMPLE:

눩 Re	sults - [Enter	. Worklist for:	DIFCF. Status: Al	I. 1 orders	Edit Moo	de:1]		0	1	100	and the second		
F 📗	ile Edit V	iew Orders	Specimens Res	ults Diffp	ad Workli	st Tool	s V	Vind	ow Help	)			
😅 🔻	· 🔒 💽 🎒	X 🖻 🖬 🚺		<b>i</b> + + + +	*   🔈 🚳	🕴 🔀 Dx	2	<b>\?</b>	🗅 🔫 🕻	💐 🕫 🎯 👘 🕹	F H		
#	TAT Status	E Pr. 🛆 Or	der Name (las	t, first, mid	dle) Tests	C Time	С	F	Patient Info	·			
<b>B</b> 1	¢	R CF	0010 TEST, CSF		CSF	15:00	2	Nε	ume (last fi	rst middle): TEST, C	SF		
				Diffpad - Template: DIFCF					Order Info	×00000059 Sex 		Wa	ard: DU№
	-		CFPOL					🔁 Refresh 🕒 Comment 😰 Coll/Rec 🛋 Verify 🛋 Verify					
			(1)	0 [2]	0 [3]				T Ord	T Name	Result	FI	Flags
								1	CSF	Appearance	Hazy		
			Total: 0					2	CSF	Color	Xanthochromic		
			Limit: 0 By:					3	CSF	CSFWBC	1236.	۱L	Н
			Ver: Y					4	CSF	CSF RBC			
			SAVED					5	CSF	Number of Cells Cour	10	0	
			SAVED										