

NWT LIS TIP OF THE WEEK

Module: SoftLabMic

Software Version: 4.0.8

Topic: Reprint Unit Conf Label

Issued On: 27 July 2018

Distribution: All Users

Page: 1 of 1

Issued By: Myriam-Laurie Dufresne, Territorial LIS Administrator

TOPIC/QUESTION:

I accidentally skipped printing my blood bank unit-retype confirmation label. How can I resend it to print?

ANSWER/TIP:

Inventory>Orders>Display

Type in Unit#...Enter>F12>Highlight the Retype Test>Enter>Ctrl-P

SOFTBANK - INVENTORY Test order

SOFTBANK - INVENTORY Test Order - Display

SOFTBANK - Unit Selection Criteria

Unit#: C052118356841 Prod: RCG (RC) ABO: O NEG Volume: 291 ml

Org. unit#: C052118356841 Src code: CBS Coll. Ctr.: C0521 Ant./attrib.: No

Label OK: Location: Comm: No Donation: allOgenic

Received: 18-07-25 16:36 by EJB2 Created/collected: 18-07-20 23:59

Status: 18-07-25 17:24 Available Expires: 18-08-31 23:59

| Order#/Test name | Description/Result | Status | Tech | Comm |
|---------------------|----------------------|-----------|------|------|
| Order#: 9125001I | req.:18-07-25 by SCC | completed | | |
| Tests: Retype Units | O NEG | completed | EJB2 | |

SOFTBANK - Unit Selection Criteria

SOFTBANK (Active)

SOFTBANK (Active) - Results of test Retype Units (RETP) view

Control

Functions

Ctrl+Q: tests

Ctrl+P: tests

F7: Logic table

Printer Name Type

View_file Y

Save_file F

R01_Lab_A_LabelPM41a L

R02_Pnl_B_LabelPM41b L

R03_Lab_PrinterL04015 P

R_03_LAND_LabelPrinter P

R11_ZEBRA_BB L

R07_ZEBRA_BB L

R04_LabBacte_PM41 L

R04_LabOutr_PM43 L

R05_LabMain_PM42 L

R06_LabYDCC_PM43 L

R04_ZEBRA_BB L

Modem

R07_Zebra L

R23_87W_Chemistry L

Unit Confirmation Label