

# NWT LIS TIP OF THE WEEK

**Module:** SoftLabMic

**Software Version:** 4.0.8

**Topic:** Glucometer Glucose  
check

**Issued On:** 9JAN2019

**Distribution:** All Users

**Page:** 1 of 1

**Issued By:** April Darrach, Territorial LIS Administrator

---

## **TOPIC/QUESTION:**

What is this add-on test that happens when I order GGTT or GTT?

## **ANSWER/TIP :**

**Test code :** GLUGL **Specimen tube type :** GLG - label will print

**Purpose:** To record the glucometer results for glucose checking prior to the administration of the Trutol drink for GGTT and GTT testing

**How to result :** It is set to mandatory at receiving and will prompt for a result as soon as specimen receiving is complete

**How it works :** the GLUGL test code will automatically reflex as soon as GGTT and GTT tests are saved on an order. Users will have to go back into Order Entry to collect and receive the GLG specimen. Glucometer results >7.8 will force another rule to cancel the GGTT and GTT test codes and will reorder a GLUF test.

For samples that are collected outside the laboratory where a glucometer check cannot be performed – enter a dot result when prompted for the glucometer result. A cancellation reason will need to be entered as well.

See **LISXXX : Glucometer Glucose Screening for GGTT and GTT** for further procedural instructions.