

# NWT LIS TIP OF THE WEEK

Module: SoftLabBank	Software Version: 4.0.8	
Topic: Clearing a PID error in Bank	Issued On:21JAN019	
Distribution: All Users	Page: 1 of 1	
Issued By: M-L Dufresne, Territorial LIS Administrators		

#### **TOPIC/QUESTION:**

I cannot access a patient record in blood bank. System just keeps telling me that the patient data is in use by another user no matter what I try to do. Where do I go from here?

### **ANSWER/TIP:**

The LIS system will enter into a locking process to prevent a second user from accessing a record when someone is already in it.

Occasionally the process gets stuck and the record remains unavailable even after the first user has exited from the record.

When this occurs:

## \*\*SCC must be engaged to kill the process\*\*

REGULAR BUSINESS HOURS		OFF HOURS	
Action	Chk	Action	Chk
Note the Process Identification Number (PID)		Note the Process Identification Number (PID)	
Contact your LIS Administrator		Phone 1-800-SoftLab and ask that process for pid= be killed for client ID:NWT	
		Advise the LIS Administrator of your call with SCC	

### **EXAMPLE:**

