

# NWT LIS TIP OF THE WEEK



<b>Module:</b> SoftLabMic/SoftBank	<b>Software Version:</b> 4.0.8
<b>Topic:</b> Sending Orders to HIS	<b>Issued On:</b> 13 February 2019
<b>Distribution:</b> All Users	<b>Page:</b> 1 of 1
<b>Issued By:</b> April Darrach, Territorial LIS Administrator	

## TOPIC/QUESTION:

When do I send orders to the HIS? How do I send orders to the HIS?

## ANSWER/TIP:

Once results are verified (Core Lab/Blood Bank) or finalized (Microbiology) the order is essentially dead so any changes made after the fact are not updated in the HIS systems. Any orders that have changes made after results have been released must have the order sent to the HIS for the changes to apply in the Viewer/EMR. To get the get the changes in the HIS' we need to "Send the order to the HIS".

Core Lab	Micro Lab	Blood Bank
<p>Perform in Order Entry</p>	<p>Perform in Micro Result Entry</p>	<p>Perform in Patient Order Screen</p>