

NWT LIS TIP OF THE WEEK

Module: SoftLabMic

Software Version: 4.0.8

Topic: Honeywell1450g Scanner
Configuration - General

Issued On: 11MAR2019

Distribution: All Users

Page: 1 of 1

Issued By: April Darrach, Territorial LIS Administrator

TOPIC/QUESTION:

How do I configure my Honeywell 1450g scanner to read the barcodes with a carriage return? This is NOT for a Blood Bank scanner

ANSWER/TIP: Scan the barcode below

WARNING: DO NOT USE THIS ON BLOODBANK SCANNERS AS THESE MUST BE CONFIGURED SPECIALLY

Barcode Symbology and Carriage Return (includes resetting the defaults):

