



Stanton Territorial Hospital
 P.O. Box 10, 550 Byrne Road
 YELLOWKNIFE NT X1A 2N1

Document Name:
 Maintenance Record – GeneXpert -1

Distribution:
 Microbiology Instrumentation Manual

Document Number:
 MIC70711

Date Issued:
DRAFT

Maintenance Record – GeneXpert-1

Month: _____

Year: _____

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
---	---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

Daily Maintenance

Disinfect work area

Weekly Maintenance

Reboot GX instrument

Reboot software and computer

Monthly Maintenance

Disinfect GX surfaces

Disinfect cartridge bay interior

Disinfect plunger rod

Clean fan filter

Archive runs

Yearly Maintenance:

Calibration of GX – performed by Cepheid

As-Required

Print system log report

Initials

Corrective action:

Reviewed by:

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME:

PRINT DATE: