

NWT LIS TIP OF THE WEEK

Module: SoftLab/SoftMic	Software Version: 4.0.8
Topic: Rejecting UCHL - Chlamydia and Gonorrhoeae, random urine	Issued On: 30 Oct 2020
Distribution: All Users	Page: 1 of 1
Issued By: Lena Schofield, Territorial LIS Administrator	

TOPIC/QUESTION:

Due to a global shortage of Aptima Urine G/C collection kits, test code: UCHL for Chlamydia and Gonorrhoeae, random urine tests need to be rejected until new stock is received. When a requisition is received ordering a UCHL, a test needs to be rejected (REJ).

ANSWER/TIP:

Order Entry

1. Complete Order Entry for the requisition. In the Ordering Grid – order **REJ** test code, **Save** Order.
2. Enter **Tests Cancelled:** UCHL
3. Enter **Reason:** P – See Below
4. Click the **Comment** button
5. Click the **F5 key** or **CanMesg (F5)** button
6. Select **@UCHL**
7. Click **OK**
8. Click **OK**
9. Click **OK**
10. A message box asking **Do you want to save the changes you made to order?** will pop-up, click **Yes**.

See the next page for pictures of the steps that need to be taken.

NWT LIS TIP OF THE WEEK

1 Select the test in the Ordered list.

2 Enter the value in the Please enter required information dialog.

3 Select the reason from the dropdown menu.

4 Click CanMesg in the Please enter required information dialog.

5 Enter the comment in the Result Comments (UNAC2)* window.

6 Select the test from the Search Results list.

7 Click OK in the Search Results dialog.

8 Click OK in the Result Comments (UNAC2)* window.

Cancelled Tests

Key	Text
A	Clotted Spec
B	Frozen Specimen
C	Hemolyzed Spec
D	Icteric Spec
E	Improper Sample
F	Insufficient Quantity
G	Leaked in Trans
H	Lipemic Spec
I	Mislabelled
J	No Sample Rec'd
K	Not Centrifuged
L	Sample too old
M	Thawed Specimen
N	Unlabelled Spec
O	report
P	See Below
Q	No smear received
R	Poor quality smear received
S	Duplicate Order

Search Results

#	ID	Text
9	@72HR	Specimen received >72 hours old
10	@7DAY	Specimen received >7 days old
11	@BF	Specimen unsuitable for analysis - fluid too viscous for
12	@CFI	Sample for INR testing was received frozen. This is NOT
13	@CGFT	Green top tube collected, unable to perform FT4
14	@DUP	Tests indicated have been ordered by multiple practition...
15	@LBL	SPECIMENS MUST BE LABELLED WITH PATIENT'S...
16	@PAN	Due to the Pandemic and as part of the Territorial Lab S...
17	@TNP	Troponin not processed on routine outpatient
18	@Uchl	APTIMA Urine Specimen Collection Kits used for Chlam...
19	@UNL	Specimen received unlabelled

Total 19 Selected 1