



Stanton Territorial Hospital
 P.O. Box 10, 550 Byrne Road
 YELLOWKNIFE NT X1A 2N1

Document Name:
 Maintenance Record – GeneXpert

Distribution:
 Microbiology PCR Testing Manual

Document Number:
 MIC81110

Date Issued:
 22 Apr 2020

Maintenance Record – GeneXpert

Month: _____ **Year:** _____

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
---	---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

Daily Maintenance

Disinfect work area

Check printer paper

Weekly Maintenance

Reboot software and laptop

Reboot GX instrument

Monthly Maintenance

Disinfect GX surfaces

Disinfect cartridge bay interior

Disinfect plunger rod

Clean fan filter(s)

Archive tests

Database management

Yearly Maintenance:

Calibration of GX – performed by Cepheid

Corrective action:

Reviewed by:

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.