

NWT LIS TIP OF THE WEEK

Module: SoftLabMic

Software Version: 4.0.8

Topic: Cancelling Type and Screen, Referrals

Issued On: 19NOV2021

Distribution: All End Users

Page: 1 of 1

Issued By: Kaylan Symes, Territorial LIS Administrator






TOPIC/QUESTION: How to properly cancel *Type and Screen, Referral* for the results to be posted to HealthNet Viewer and Electronic Medical Records.

ANSWER/TIP : Cancelling of bloodbank orders is to be performed by an MLT and *only in the SoftBank Module*.

Always cancel TSR by the group test, not the individual tests as this will not flow to EMR/HNV

EXAMPLE :

1. Open SoftBank
2. Click 1. Patient
3. O Orders
4. New section opens, click Cancel
5. Enter patient information, press F12
6. Select correct stay, press enter
7. F12 to accept the patient information
8. Highlight the *group test*, press enter
9. Pop up appears, click yes
10. Enter cancellation comments, click F12



Order will appear in SoftLab and SoftBank :

T Ord	T Ind	Result	FI	Flags	M	Prev Res	Prev Date	C	Status	Tech	DateV	TimeV	TechV
1	TSR	ABOR							LC V	KDS	2021-11-17	15:35	KDS
2	TSR	RHR							LC V	KDS	2021-11-17	15:35	KDS
3	TSR	ABSCR							LC V	KDS	2021-11-17	15:35	KDS

Type	ID	Priority	Cycl
Ordered (0)			
Cancelled (1)			

Order/Test/Product	Description/Result	Status	Tech	Comm
Order#: G7170004	req.:21-11-17 by TESD1			
Tests: TS Referral		cancelled		!C
ABO Referral		cancelled	KDS	
Rh Referral		cancelled	KDS	
Ab Screen Ref		cancelled	KDS	