
 <p><b>Stanton Territorial Hospital</b> P.O. Box 10, 550 Byrne Road YELLOWKNIFE NT X1A 2N1</p>	<p><b>Document Name:</b> Labelling Blood Culture Bottles Job Aid</p> <p><b>Distribution:</b> Microbiology Culture Manual</p>	<p><b>Document Number:</b> MIC35100</p> <p><b>Date Issued:</b> November 8, 2021</p>
---	--	---

1. Ensure the correct specimen label is placed on the correct blood culture bottle:

	<p><b>BAE</b> label is placed on <b>Aerobic Bottle</b></p>
	<p><b>BAN</b> label is placed on <b>Anaerobic Bottle</b></p>
	<p><b>BPE</b> label is placed on <b>Pediatric Bottle</b></p>

2. Ensure the specimen label is placed on the blood culture bottle correctly. Do not place LIS specimen label ovetop of the bottle barcode:



**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.