



Stanton Territorial Hospital
P.O. Box 10, 550 Byrne Road
YELLOWKNIFE NT X1A 2N1

Document Name:
Maintenance Record – Vitek 2

Distribution:
Microbiology Instrumentation Manual

Document Number:
MIC70310

Date Issued:
17 Nov 2017

Maintenance Record – Vitek 2

Month: _____ **Year:** _____

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
---	---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

Daily Maintenance

Empty waste tray																															
Check status																															
Zero DensiCHEK																															
Check printer paper																															
Initials																															

Monthly Maintenance

Task	Date	Tech	Task	Date	Tech	Task	Date	Tech
Clean boats			Clean vacuum seal			Optics diagnostic test		
Clean carousels			Clean vacuum chamber			Clean Vitek 2 exterior		
Clean cassettes			Clean optics			Clean Vitek 2 computer		
Clean waste tray			Change saline bag			Clean Smart Carrier		
Clean base pan			Change pipette tips			Clean DensiCHEK		
Clean drip pan			Disp/pipette volume test			Sterilize dispensette		

DensiCHEK plus Calibration

Calibration Standard	Reference Range	Lot #	Expiry	Result	Date	Tech
0.5	0.44 to 0.56					
2.0	1.85 to 2.15					
3.0	2.79 to 3.21					

Corrective action: _____ Reviewed by: _____

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.