



**Stanton Territorial Hospital**  
 P.O. Box 10, 550 Byrne Road  
 YELLOWKNIFE NT X1A 2N1

**Document Name:**  
 Maintenance Record – GeneXpert

**Distribution:**  
 Microbiology PCR Testing Manual

**Document Number:**  
 MIC81110

**Date Issued:**  
 22 Apr 2020

**Maintenance Record – GeneXpert** Instrument: \_\_\_\_\_ Month: \_\_\_\_\_ Year: \_\_\_\_\_

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
---	---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

**Daily Maintenance**

Disinfect work area

Check printer paper

**Weekly Maintenance**

Reboot software and laptop

Reboot GX instrument

**Monthly Maintenance**

Disinfect GX surfaces

Disinfect cartridge bay interior

Disinfect plunger rod

Clean fan filter(s)

Archive tests

Database management

**Yearly Maintenance:**

Calibration of GX – performed by Cepheid

Corrective action:

Reviewed by:

\_\_\_\_\_

**NOTE:** This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC81110-Maintenance Record-GeneXpert.docxFOR.doc

PRINT DATE: 5/28/2020 9:40:00 AM