



Stanton Territorial Hospital
 P.O. Box 10, 550 Byrne Road
 YELLOWKNIFE NT X1A 2N1

Document Name:
 Maintenance Record-GeneXpert

Distribution:
 Microbiology PCR Testing Manual

Document Number:
 MIC81110

Date Issued:
 22 Apr 2020

Maintenance Record – GeneXpert Instrument: _____ Month: _____ Year: _____

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
---	---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

Daily Maintenance

Disinfect work area																															
Check printer paper																															

Weekly Maintenance

Reboot software and laptop																														
Reboot GX instrument																														

Monthly Maintenance

Disinfect GX surfaces																														
Disinfect cartridge bay interior																														
Disinfect plunger rod																														
Clean fan filter(s)																														
Archive tests																														
Database management																														

Yearly Maintenance:

Calibration of GX – performed by Cepheid																														
--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

Corrective action:

Reviewed by:

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.