



Uncontrolled When Printed

1. Ensure the correct specimen label is placed on the correct blood culture bottle:

	<p>BAE label is placed on Aerobic Bottle</p>
	<p>BAN label is placed on Anaerobic Bottle</p>
	<p>BPE label is placed on Pediatric Bottle</p>

2. Ensure the specimen label is placed on the blood culture bottle correctly. Do not place LIS specimen label ovetop of the bottle barcode:



NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.