

Change black air filter

Stanton Territorial Hospital

P.O. Box 10, 550 Byrne Road YELLOWKNIFE NT X1A 2N1

Document Name:

Maintenance Record-BACTEC FX

Distribution:

Microbiology Instrumentation Manual

Document Number:

MIC711110

Date Issued:

10 May 2019

Clean exterior surfaces

Reviewed by:

Maintenance Record - BACTEC FX						-	Month:								Year:																
	1	2	3	4	5	6	7	8	9	1 0	1	1 2	1 3	1 4	1 5	1 6	1 7	1 8	1 9	2	2	2 2	2	2 4	2 5	2	2 7	2	2 9	3	3 1
Daily Maintenance - Left	Side	e (A	1)																												
Digital Temp: 33.5 to 36.5																															
Internal Temp: 33.5-36.5																															
Green LEDs																															
Red LEDs																															
Yellow LEDs																															
Daily Maintenance - Rigl	nt Si	de ((B)																												
Digital Temp: 33.5-36.5																															
Internal Temp: 33.5-36.5																															
Green LEDs																															
Red LEDs																															
Yellow LEDs																															
External																															
Green lights																															
Red lights																															
Yellow lights																															
Audible alarm																															
Check current inventory																															
Initials																															
Monthly Maintenance																															
Task	Dat	te		Ted	ch			Ta	ask				D	ate		T	ech				Tas	sk				D	ate		-	Tecl	h

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not co	ontrolled and should be checked against electronic version prior to
use.	ontrolled and enough by encoured against electrolle version prior to
use.	
FILENAME: MIC71110-Maintenance Record-BACTEC FX	PRINT DATE: 9/30/2024 8:35:00 AM

Change silver air filter

Corrective action: