NORTHWEST TERRITORIES	Stanton Territorial Hospital	Document Name:	Document Number:
	P.O. Box 10, 550 Byrne Road	Maintenance Record-GeneXpert System	MIC72110
Health and Social	YELLOWKNIFE NT X1A 2N1	Distribution:	Date Issued:
Services Authority		Microbiology Instrumentation Manual	22 Apr 2020

<u> Maintenance Record – GeneXpert System</u>								Instrument:						Month:								Year:												
	1		2	2	3	4	5	6	7	8	9	1 0	1 1	1 2	1 3	1 4	1 5	1	1	L 7	1 8	1 9	2 0	2 1	2 2	2 3	2 4	2 5	2 6	2 7	2	2 8	2 9	33 01
Daily Maintenance																																		
Disinfect work area																																		
Check printer paper																																		
Initials																																		
						1																-										-		
	1	1	2	3	4	5	6	7	7 8	3 9	9	1 : D :	1 1	1 2	1 3	1 4	1 5	1 6	1 7	1 8	1 9	2	1	2	2 3			5	6	2 7	2 8	2 9	3 0	3 1
Weekly Maintenance																																		
Reboot software																																		
Reboot laptop																																		
Reboot GX instrument																																		
Initials																																		

Monthly Maintenance														
Task	Date	Tech	Task	Date	Tech	Task	Date	Tech						
Disinfect GX surfaces			Disinfect plunger rods			Archive tests								
Disinfect cartridge bay			Clean fan filter(s)			Database management								

Corrective action:

Reviewed by:

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC72110-Maintenance Record-GeneXpert System System

PRINT DATE: 10/1/2024 10:58 AM