



Stanton Territorial Hospital
 P.O. Box 10, 550 Byrne Road
 YELLOWKNIFE NT X1A 2N1

Document Name:
 Maintenance Record-BIOFIRE TORCH

Distribution:
 Microbiology Instrumentation Manual

Document Number:
 MIC73110

Date Issued:
 22 Apr 2020

Maintenance Record – BIOFIRE TORCH

Month:

Year:

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
---	---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

Daily Maintenance

Disinfect work area																															
Check printer paper																															
Initials																															

Monthly Maintenance

Task	Date	Tech	Task	Date	Tech	Task	Date	Tech
Disinfect Torch top and sides			Disinfect Torch screen			Clean used fan filters		
Disinfect Torch module covers			Change fan filters			Archive tests		

Corrective action: _____ Reviewed by: _____

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.