

 <b>Stanton Territorial Hospital</b> P.O. Box 10, 550 Byrne Road YELLOWKNIFE NT X1A 2N1	<b>Document Name:</b> Labelling Blood Culture Bottles Job Aid  <b>Distribution:</b> Microbiology Specimen Processing Manual	<b>Document Number:</b> MIC10180  <b>Date Issued:</b> November 8, 2021
--	---	--

**Uncontrolled When Printed**

1. Ensure the correct specimen label is placed on the correct blood culture bottle:

	<b>BAE</b> label is placed on <b>Aerobic Bottle</b>
	<b>BAN</b> label is placed on <b>Anaerobic Bottle</b>
	<b>BPE</b> label is placed on <b>Pediatric Bottle</b>

2. Ensure the specimen label is placed on the blood culture bottle correctly. Do not place LIS specimen label ovetop of the bottle barcode:



**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC10180-Labeling Blood Culture Bottles Job Aid

Print Date: 11/7/2021 2:35 PM