



**Stanton Territorial Hospital**  
P.O. Box 10, 550 Byrne Road  
YELLOWKNIFE NT X1A 2N1

**Document Name:**  
Maintenance Record-BACTEC FX

**Distribution:**  
Microbiology Instrumentation Manual

**Document Number:**  
MIC711110

**Date Issued:**  
10 May 2019

**Maintenance Record – BACTEC FX**

**Month:**

**Year:**

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
<b>Daily Maintenance – Left Side (A)</b>																															
Digital Temp: 33.5 to 36.5																															
Internal Temp: 33.5-36.5																															
Green LEDs																															
Red LEDs																															
Yellow LEDs																															
<b>Daily Maintenance – Right Side (B)</b>																															
Digital Temp: 33.5-36.5																															
Internal Temp: 33.5-36.5																															
Green LEDs																															
Red LEDs																															
Yellow LEDs																															
<b>External</b>																															
Green lights																															
Red lights																															
Yellow lights																															
Audible alarm																															
Check current inventory																															
<b>Initials</b>																															

<b>Monthly Maintenance</b>								
<b>Task</b>	<b>Date</b>	<b>Tech</b>	<b>Task</b>	<b>Date</b>	<b>Tech</b>	<b>Task</b>	<b>Date</b>	<b>Tech</b>
Change black air filter			Change silver air filter			Clean exterior surfaces		

Corrective action: \_\_\_\_\_ Reviewed by: \_\_\_\_\_

**NOTE:** This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC711110-Maintenance Record-BACTEC FX

PRINT DATE: 5/15/2025 7:43:00 PM