



## TRAINING UPDATE

**Lab Location:** SGAH & WAH  
**Department:** Core

**Date Distributed:** 3/20/2013  
**Due Date:** 4/1/2013

### DESCRIPTION OF PROCEDURE REVISION

<b>Name of procedure:</b>
<b>Advanced Micro-Osmometer Model 3320 Analyzer Maintenance Log</b> <b>Form # AG.F161.001</b>
<b>Description of change(s):</b>
Added quarterly calibration as reminder of requirement by the SOP (form AG.F110 is still required to document actual calibration)  Moved Air Vent and Solenoid to 'as needed maintenance'  Changes are shown in yellow highlight on attached form. Revised form will be put into use on 4/1/13.

Document your compliance with this training update by taking the quiz in the MTS system.

