

**TRAINING UPDATE**

**Lab Location:** SGAH and WAH      **Date Implemented:** 12.19.2013  
**Department:** Blood Bank      **Due Date:** 12.31.2013

**DESCRIPTION OF PROCEDURE REVISION**

<b>Name of procedure:</b>
Daily Temperature Quality Control
<b>Description of change(s):</b>
<ul style="list-style-type: none"><li>• Removed undercounter refrigerator</li><li>• Added space to document monthly maintenance of plasma thawer</li></ul>



### Daily Temperature Quality Control

<b>Month:</b>	<b>Year:</b>	<b>Reviewed:</b>
---------------	--------------	------------------

Date	Heat Block SN 406N0367 (36 -38C)	Hole	Heat Block SN 410N066 (36-38C)	Hole	Plasma Thawer SN 51351-802 (30-37C)	Room Temp (20-24C)	Interp S or N	Tech
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
13								
14								
15								
16								
17								
18								
19								
20								
21								
22								
23								
24								
25								
26								
27								
28								
29								
30								
31								

Date	Tech	
		<b>Plasma Thawer Water Changed</b>

**Comments/Corrective Action:**

Document: AG.F236[1] Status: RELEASED, Effective: 12/19/2013, Check Version Before Use