

TRAINING UPDATE

Lab Location: SGAH & WAH
Department: Micro

Date Distributed: 11/13/2014
Due Date: 12/1/2014
Implementation: **12/1/2014**

DESCRIPTION

Name of document:
Room Temperature Log, Microbiology AG.F314.0
Description of change(s):
<p>New Room Temperature log specifically for Micro bench / area.</p> <ul style="list-style-type: none">• The acceptable range is DIFFERENT from what is used in all other areas of the lab• A min / max thermometer is required <p>This FORM will be implemented on December 1, 2014</p>

Document your compliance with this training update by taking the quiz in the MTS system.



- Shady Grove Adventist Hospital
- Washington Adventist Hospital
- Germantown Emergency Center

Room Temperature Log - Microbiology

Acceptable range: 15 - 30C

Month: _____	Year: _____	Thermometer ID: _____
--------------	-------------	-----------------------

Interp: S = satisfactory, temp w/in acceptable range
 U = unsatisfactory, circle result and document corrective action; notify supervisor/Incharge immediately.
 Complete Variance form. Additional corrective action should be documented on the back of the sheet.
 Check Reset box to document device is reset daily

Date	Minimum Temp.	Maximum Temp.	Interp	Reset	Tech Code	Corrective action	Weekly Review
1							
2							
3							
4							
5							
6							
7							
8							
9							
10							
11							
12							
13							
14							
15							
16							
17							
18							
19							
20							
21							
22							
23							
24							
25							
26							
27							
28							
29							
30							
31							

Monthly review: _____

Date: _____

Corrective Action Codes:

- | | |
|---|---|
| <p>A. Rechecked temperature <i>(Log repeat temperature)</i></p> <p>B. Adjusted Thermostat <i>(Retake temperature in 1 hour)</i></p> <p>C. Contacted Plant Ops / Facilities</p> | <p>D. Affected contents relocated</p> <p>E. Contents quarantined pending evaluation</p> <p>F. Other: <i>(Document on reverse side)</i></p> |
|---|---|