

TRAINING UPDATE

Lab Location: SGAH and WAH **Date Implemented:** 12.22.2014
Department: Blood Bank **Due Date:** 01.15.2015

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:

Manual Product Storage Temperature Form

Description of change(s):

1. Added reason for starting manual temps to form.
2. Changed GEC acceptable range to 1-6°C.
 - a. RhIG is no longer stored in the GEC lab and all of their reagents may be stored at this temperature range.
 - b. BB will remain 2-6°C due to our reagent requirements.

Electronic Document Control System



Document No.: AG.F51[1]

Title: BLOOD BANK MANUAL PRODUCT STORAGE TEMPERATURE FORM

Owner: LESLIE.X.BARRETT LESLIE BARRETT

Status INWORKS

Effective Date: 21-Jan-2015

Next Review Date:

