

## TRAINING UPDATE

**Lab Location:** GEC  
**Department:** Core

**Date Distributed:** 6/2/2015  
**Due Date:** 6/30/2015  
**Implementation:** 6/2/2015

### DESCRIPTION OF PROCEDURE

<b>Name of procedure:</b>
<b>Dimension Xpand Daily System Check    AG.F290.1</b>
<b>Description of change(s):</b>
<p>Replace ABS product with CHK and update appropriate values</p> <p><b>This revised FORM was implemented on June 2, 2015</b></p>

Document your compliance with this training update by taking the quiz in the MTS system.

**DIMENSION® XPAND DAILY SYSTEM CHECK**

Instrument Serial #: \_\_\_\_\_

Month &amp; Year \_\_\_\_\_

All wavelengths OK?	Reagent 1		Reagent 2		Sampler		HM		CHK Lot	Tech
	Mean	SD	Mean	SD	Mean	SD	Mean	SD		
293 filter: -2.5 to +2.5 All other filters: -1.5 to +1.5	CHK Carton Value ± 15	≤ 3.8	CHK Carton Value ± 15	≤ 3.8	10% of CHK Carton Value ± 2	≤ 1.6	10% of CHK Carton Value ± 4	≤ 1.6		
1										
2										
3										
4										
5										
6										
7										
8										
9										
10										
11										
12										
13										
14										
15										
16										
17										
18										
19										
20										
21										
22										
23										
24										
25										
26										
27										
28										
29										
30										
31										

Weekly review:	Weekly review:	Weekly review:
Weekly review:	Weekly review:	Monthly review: