

TRAINING UPDATE

Lab Location: GEC, SGMC & WAH
Department: Core

Date Distributed: 1/5/2016
Due Date: 1/31/2016
Implementation: 2/1/2016

DESCRIPTION OF PROCEDURE REVISION

Name of procedure:
Correcting Results GEC.LIS45, SGAH.LIS52, WAH.LIS51 v1
Description of change(s):
<p>Section 4: add UCR, CLR, CREAT Section 5:</p> <ul style="list-style-type: none">• add criteria for corrections (part A)• specify process (must press <i>enter</i> through each result to assure calculations will also correct) and add examples to step B.2• edit format for call documentation (same as for documenting a critical value call)• add note for correcting panels (must add comment to EVERY result that was changed)• replace PI form with QV <p>This revised SOP will be implemented on February 1, 2016</p>

Document your compliance with this training update by taking the quiz in the MTS system.

Approved draft for training (version 1)

Non-Technical SOP

Title	Correcting Results	
Prepared by	Marie Sabonis	Date: 12/21/2009
Owner	Marie Sabonis	Date: 12/21/2009

Laboratory Approval		
Print Name and Title	Signature	Date
<i>Refer to the electronic signature page for approval and approval dates.</i>		
Local Issue Date:		Local Effective Date:

Review:		
Print Name	Signature	Date

TABLE OF CONTENTS

1. PURPOSE.....	2
2. SCOPE	2
3. RESPONSIBILITY.....	2
4. DEFINITIONS.....	2
5. PROCEDURE.....	3
6. RELATED DOCUMENTS	7
7. REFERENCES	7
8. REVISION HISTORY.....	8
9. ADDENDA AND APPENDICES.....	8

1. PURPOSE

Whenever an error on accepted (released) results is discovered, the corrected results are called to the appropriate area. This procedure describes the steps to correct a result in the LIS and enter the call documentation.

2. SCOPE

Any reports that require correction will be processed via this procedure. No results should ever be deleted. The correction process documents the previous report as a comment after the corrected result.

3. RESPONSIBILITY

All technical staff will maintain current knowledge of this function.

4. DEFINITIONS

OEM	Online Result Entry
MEM	Manual Result Entry
LIS	Laboratory Information System
M-ACC#	Modify Accession Number Prompt
WBC	White Blood Cell
NRBC	Nucleated Red Blood Cell
A/M/D/P/R	Accept/ Modify/ Display/ Preliminary/ Reject
CRW	Credit Test Request without Update of Patient Data
UCRR	Urine Creatinine
CLR	Creatinine Clearance
CREAT	Creatinine
BMP	Basic Metabolic Profile

5. PROCEDURE

A. Criteria for Correcting Results

1. Results that should be corrected
 - Values that change from normal (within the reference range) to abnormal or vice versa
 - Difference between values exceeds the TEA
 - Results do not belong to the patient
2. Results that do NOT need correction
 - Both values are within reference range (normal)
 - Both values are outside the reference range (abnormal)
 - Difference between values is less than TEA

B. Core Laboratory

This procedure is to be used to correct any reported erroneous result. DO NOT CREDIT OR CHANGE REPORTS VIA LIS FUNCTION OER. Using either function OEM, keyboard, or keyboard MEM, the Medical Technologist/Technician corrects the result. The individual correcting the result will take the following steps:

1. **Call the appropriate nursing unit, outpatient area, emergency department, or physician, as soon as the error has been discovered, and report the error to the nurse, physician assistant, or physician caring for the patient.**
2. Access the accession number in function MEM or OEM, using M-ACC# to correct the resulted test.
Note: When modifying ANY result(s) you MUST press *ENTER* through each result field until you reach the 'ACCEPT, MODIFY, REJECT?' prompt before making any changes. Once at the prompt, type M to Modify and make the appropriate changes.

This process is required to assure that any calculations associated with the corrected test results will be re-calculated. See examples below and refer to step 4 for proper call documentation.

a. Example of correct process

The UCRR result needed to be changed.

ENTER was pressed until the ACCEPT, MODIFY, REJECT prompt was reached.

M (Modify) was typed and *ENTER* pressed until test result that needed to be changed.

The correct value was keyed and *ENTER* pressed through the remaining prompts. Note that the calculated test CLR has been re-calculated.

DEVICE	LOC:	WAH	WASHINGTON ADVENTIST HOSPITAL	HOSP. ID:	WAH	
ACC NO	NAME	PN:	TEST-8	AGE/SEX	LOC	PHYSICIAN
T689	TEST, MARIE			35Y F	TEST	CACCIABEVE MD,
			DOB: 10/02/1980		COLL: 10/20/2015 12:18	

CREAT : 0.2,*,L
 UCRR : 1.5,*,L
 2.0 FAILED NORMAL [30.0-125.0]mg/dL
 CORRECTED ON 10/20 AT 1303: PREVIOUSLY REPORTED AS 1.5
 HT : 63
 WT3 : 125
 CLR : 7,*,L
 9
 **CALC WAS REAPPLIED USING MODIFIED RESULTS FOR UCRR **
 FAILED NORMAL [80-120]mL/min/m2
 CORRECTED ON 10/20 AT 1303: PREVIOUSLY REPORTED AS 7
 PER : 24

b. Example of *incorrect process*

The BMP was resulted and the CREAT result needed to be changed. After typing in M-ACC#, *ENTER* was pressed until the CREAT test was reached and the CREAT result was changed. Then *ENTER* was pressed until the ACCEPT, MODIFY, REJECT prompt was reached. Note that the GFR calculation did **NOT** update as it should have since the creatinine value is used for the calculation of the GFR.

MANUAL RESULT ENTRY

ACC. NO.: M-T691

SOD: 147,*,H

K: 4.2

CL: 102

CO2: 23

BUN: 47,*,H

CREAT: 0.2,*,L
 0.4 FAILED NORMAL [0.6-1.3]mg/dL
 FAILED DELTA: 0.2 10/20/2015 1218 T689
 ACCEPT (Y/N)? : Y
 CORRECTED ON 10/20 AT 1351: PREVIOUSLY REPORTED AS 0.2

GLUC: 102

CA: 8.9

AGAP: 26,*,H
 26 FAILED NORMAL [5.0-15.0]

GFR: >90
 NO NEW RESULTS FOUND FOR GFR
 UNABLE TO CALCULATE GFR, NO CHANGES FILED FOR GFR

GFRA: >90
 NO NEW RESULTS FOUND FOR GFRA
 UNABLE TO CALCULATE GFRA, NO CHANGES FILED FOR GFRA

(ORDER: BMP)

Form revised 3/31/00

ACCEPT (A), MODIFY (M), DISPLAY PRIOR (D), PRELIM (P), OR REJECT (R)?

Note: The ACCEPT, MODIFY, REJECT prompt **MUST** be reached before selecting M (Modify) to make changes. This ensures that any applicable calculations will update.

3. The laboratory corrects results for two primary reasons:
 - (1) to replace the incorrect results with the correct results, a simple exchange of results, or
 - (2) when we have been notified that the tested specimen did not belong to the patient. In this situation, the incorrect results are replaced with an English Text code such as WPT (Wrong patient tested, please disregard these results).
4. For both corrections, the computer will generate an automatic message, "Previously reported as xxx" which leaves the incorrect results attached to the patient's chart.

a. To enter correct numeric results:

- 1) With cursor flashing underneath the incorrect result, type in the correct result, and append the call documentation as **free text with the full name of the person taking the report and the date and time.**

GLU: 100
You Type 180-; **SUE SMITH 110515 1830**

Note: When correcting multiple results within a panel (such as BMP or CBC) call documentation is appended to each result that was changed.

- 2) After pressing the enter key, the system will append the following messages to the result:

GLU: 180
SUE SMITH 110515 1830
CORRECTED ON 01/06/02 AT 0930:
PREVIOUSLY REPORTED AS 100

- 3) **Accept** the change at the A/M/D/P/R prompt.

b. To enter English Text code in place of erroneous numeric values:

- 1) With cursor flashing in the result entry field, type in the English text code and append the documentation of who was called and the time.

GLU: 100
You Type WPT-; **SUE SMITH 110515 0930**

- 2) When enter is pressed, the system will append the following statement to the result:

GLU: *(No result will show beside the test code due to the English Text Code)*

Wrong Patient Tested, please disregard these results, SUE SMITH 110515 0930.

CORRECTED ON 01/06/02 AT 0930:

PREVIOUSLY REPORTED AS 100

- 3) **Accept** the change at the A/M/D/P/R prompt.

5. Complete a [Quality Variance \(QV\)](#) form and attach a copy of the revised report.

C. Microbiology

This procedure is to be used to correct any finalized erroneous microbiology result. The Medical Technologist/Technician will correct the results in the Microbiology Result Entry GUI application. The LIS will automatically generate a corrected report statement to attach to new result.

Note: Changes made to preliminary results will not generate a corrected statement.

Changes to quantity will not generate a corrected statement.

The individual correcting the result will take the following steps:

1. Call the appropriate nursing unit, outpatient area, emergency department, or physician, as soon as the error has been discovered, and report the error to the nurse, physician assistant, or physician caring for the patient.
2. Enter the corrected result.
 - a. From the Gateway screen select **Microbiology Result Entry** by double clicking or highlight it and press **enter**.
 - b. Enter the Accession number and press **enter**.
 - c. Verify the correct test is highlighted and press **enter** to select it.
 - d. In the Microbiology Result Entry screen select the correct tab for the test that needs to be corrected.
 - 1) Direct Exam- Wet Prep, Gram Stain, Malaria (ALT+D)
 - 2) Culture Entry- Blood Culture (ALT+L)
 - 3) Susceptibility- Drug sensitivities (ALT+U)
 - e. Use the **up-arrow** key to move the cursor over the incorrect result. The result will have a blue box around it.
 - f. Press the key for the correct result. (If you don't know the correct result key, press F8 to make the keyboard visible and select the key from there)

8. REVISION HISTORY

Version	Date	Reason for Revision	Revised By	Approved By
		Supersedes SOP LIS039.002		
000	12/10/15	Section 4: add UCR, CLR, CREAT Section 5: add criteria for corrections (part A), specify process and add examples to step B.2, edit format for call documentation, add note for correcting panels, replace PI form with QV Footer: version # leading zero's dropped due to new EDCS in use as of 10/7/13.	M Sabonis	M Sabonis

9. ADDENDA AND APPENDICES

None