

TRAINING UPDATE

Lab Location:
Department:

GEC, SGMC & WAH
Core

Date Distributed:

2/19/2016

Due Date:

2/29/2016

Implementation:

3/1/2016

DESCRIPTION OF REVISION

Name of procedure:
HCG Quality Control Log AG.F42 v2
Description of change(s):
Added column to indicate if sample is urine or serum This revised FORM will be implemented on March 1, 2016

Document your compliance with this training update by taking the quiz in the MTS system.

HCG QUALITY CONTROL LOG

- Germantown Emergency Center
- Shady Grove Medical Center
- Washington Adventist Hospital

1. **External Controls** are tested and documented once each week and with each new kit lot number or shipment.
2. Internal controls must be documented each time the test is performed.

Date	Sample Type	Patient Name / MR#	Patient Result	Kit	Internal Control	External Pos Control		External Neg Control		Tech
				Lot # / Expire	Valid / Invalid	Lot # / Expire	Result	Lot # / Expire	Result	
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
Weekly review:				Weekly review:			Weekly review:			
Weekly review:				Weekly review:			Monthly review:			