

TRAINING UPDATE

Lab Location: GEC
Department: Core Lab

Date Distributed: 7/30/2018
Due Date: 8/31/2018
Implementation: 8/1/2018

DESCRIPTION OF REVISION

Name of procedure:
i-STAT 1 System Maintenance Log AG.F213.2
Description of change(s):
<p>Removed section for 6 month calibration verification. Why? This is check is not performed by bench staff; it is performed by technical specialist and documented in the Recurring Schedule binder</p> <p>This revised form was implemented on August 1, 2018</p>

Document your compliance with this training update by taking the quiz in the MTS system.

- Germantown Emergency Center
- Shady Grove Medical Center
- Washington Adventist Hospital

i-STAT 1 System Maintenance Log

Month: _____

Year: _____

Instrument Serial Number: _____

Daily	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Check Printer Paper																															
Decontaminate the Bench and instrument exterior																															
Run the Electronic Stimulator																															
Tech																															

Weekly					
Run Blood Gas QC	Date:	Date:	Date:	Date:	Date:
	Tech:	Tech:	Tech:	Tech:	Tech:
Run BNP QC	Date:	Date:	Date:	Date:	Date:
	Tech:	Tech:	Tech:	Tech:	Tech:

As Needed	
Replace the Batteries	Date:
	Tech:

Date	Problem	Action Taken	Tech

Weekly review:	Weekly review:	Weekly review:
Weekly review:	Weekly review:	Monthly review:

Document: AG.F213[2] Status: PRERELEASED, Effective: 1/1/2099, Check Version Before Use