

## TRAINING UPDATE

**Lab Location:** SGMC & WAH  
**Department:** Core

**Date Distributed:** 3/7/2019  
**Due Date:** 3/31/2019  
**Implementation:** 4/1/2019

### DESCRIPTION OF REVISION

**Name of procedure:**

**Iris Maintenance Logs AG.F46.4**

**Description of change(s):**

Removed Body Fluid QC and background checks from page 1 (testing now done on Sysmex).

Updated logo on both pages.

**This revised FORM will be implemented on April 1, 2019**

**Document your compliance with this training update by taking the quiz in the MTS system.**

- Germantown Emergency Center
- Shady Grove Medical Center
- Washington Adventist Hospital

Month \_\_\_\_\_ Year \_\_\_\_\_

Day of the Month	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
------------------	---	---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

**EACH SHIFT**

** Run Urine Controls– day tech																																	
** Run Urine Controls – eve tech																																	
** Run Urine Controls – night tech																																	

**DAILY**

Clean Instrument Surfaces																																
Clean the Sampler																																
Clean Load/Unload Station (optional)																																
Check Lamina Supply																																

**MONTHLY**

Perform Instrument Calibration																																
Perform Backup																																

**AS NEEDED**

Clean Rinse/Waste Bath																																
Replace Lamina Container																																
Clean Sample Tube Detector																																
Clean Barcode Reader Window																																
Inspect and Clean Racks																																
Clean Optical Sensors																																
Replace Lamina Filter																																
Clean Sample Filter																																

Weekly review:	Weekly review:	Weekly review:
Weekly review:	Weekly review:	Monthly review:

\*\* Always run positive QC prior to the negative control to verify carry over has not occurred.

- Germantown Emergency Center
- Shady Grove Medical Center
- Washington Adventist Hospital

### AX-4280 Maintenance Log

Month \_\_\_\_\_ Year \_\_\_\_\_

DAILY	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Empty the Waste Box																															
Empty Liquid Waste																															
Check Strip Supply																															
Check Wash Supply																															
Run Controls																															
Tech																															

WEEKLY				
Date				
Clean SG Cell				
Clean the Test Strip Feeder				
Clean the Test Strip Tray				
Calibration Verification				
Perform SG Calibration				
Tech				

MONTHLY	
Date	
Clean Waste Filters	
Clean Washing Bath	
Tech	

QUARTERLY	
Date	
Clean/Replace Wash Filters	
Tech	

AS NEEDED	
Date	
Clean Trap Bottle	
Tech	

Weekly review:	Weekly review:	Weekly review:
Weekly review:	Weekly review:	Monthly review: