

## TRAINING UPDATE

**Lab Location:** SGMC & WAH  
**Department:** Core

**Date Distributed:** 7/1/2019  
**Due Date:** 7/31/2019  
**Implementation:** 8/1/2019

### DESCRIPTION OF REVISION

**Name of procedure:**

**Sysmex Stainer SP-10 Maintenance Log AG.F378.1**

**Description of change(s):**

- Changed layout of form from landscape to portrait
- Changed heading for tech to Tech Code
- Updated logo

There is no change to actual content of the log

**This revised FORM will be implemented on August 1, 2019**

**Document your compliance with this training update by taking the quiz in the MTS system.**

Month / Year \_\_\_\_\_

Instrument Serial Number: \_\_\_\_\_

Date	Perform Shutdown 1 √	Clean Spreader Glass √	Check Trap Chamber √	Clean Stain Cassettes √	Tech Code
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					
16					
17					
18					
19					
20					
21					
22					
23					
24					
25					
26					
27					
28					
29					
30					
31					

Weekly	Week 1	Week 2	Week 3	Week 4	Week 5
Perform Shutdown 2 √					
Tech Code					
Date					

Monthly	√	Tech Code	Date
Clean Racks			
Clean right and left sampler rack pool			
Clean Measurement Line			

Action / Comment: \_\_\_\_\_  
 \_\_\_\_\_

Weekly review:	Weekly review:	Weekly review:
Weekly review:	Weekly review:	Monthly review: