

## TRAINING UPDATE

**Lab Location:** FWMC  
**Department:** Core Lab

**Date Distributed:** 11/9/2021  
**Due Date:** 12/9/2021

### DESCRIPTION OF REVISION

<b>Name of procedure:</b>
AG.FW20.2 Sure-View Serum and Urine HCG STAT Test QC Log
<b>Description of change(s):</b>
<b>Added column to record the sample type (urine or serum)</b>
<b>This form will be implemented on November 15, 2021</b>

Document your compliance with this training update by taking the quiz in the MTS system.

1. **External Controls** are performed each day of patient testing.
2. **Internal controls** must be documented each time a patient test is performed.
3. If QC results are not acceptable, document corrective action. Do not accept patient results before reviewing QC results for proper reactions.

Date	Sample Type	Patient Name / MR#	Patient Result	Kit	Internal Control	External Positive Control		External Negative Control		Tech Code
				Lot # / Expire	Red Control Line (Y / N)	Lot # / Expire	Result	Lot # / Expire	Result	
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
	<input type="checkbox"/> Serum <input type="checkbox"/> Urine									
Weekly review:				Weekly review:			Weekly review:			
Weekly review:				Weekly review:			Monthly review:			