

Laboratories at Adventist Healthcare White Oak Medical Center, Shady Grove Medical Center, Fort Washington Medical Center and Germantown Emergency Center

Date: 2/11/22 LABORATORY ALERT Subject:

The Laboratory is informing you of the following change:

Effective date:	February 15, 2022
SOP	BioFire® FilmArray® Meningitis/Encephalitis (ME) Panel (AHC.M1022)
Test Name:	Meningitis/Encephalitis (ME) Panel
Test Code:	BFME
Specimen Requirements:	CSF that <u>has not been spun</u> . If the tube used for Chemistry testing (glucose and protein) and has been spun down – we cannot use this tube to run on the BioFire because it could cause false results.
Additional information	An MTS has been assigned to SGMC and WOMC technical and non- technical staff to make sure that both lab sections are aware that spun CSF cannot be used as an "add on" for this test. If you have any questions please ask your lead tech or supervisor.
Contacts:	Stephanie Codina, Laboratory Manager 240 6015 Hollie Genser, SGMC Lab Supervisor 240-826-6432 Zanetta Morrow, WOMC Lab Supervisor 240-637-5681 Derrick Wallace, SGMC Field Ops Supervisor 240-826- 6642