

**TRAINING UPDATE**

**Lab Location:** All Sites      **Date Implemented:** 9/26/23  
**Department:** Blood Bank      **Due Date:** 10/15/23

**DESCRIPTION OF PROCEDURE REVISION**

<b>Name of procedure:</b>
Infrared Thermometer Quality Control (Fisher Traceable)
<b>Description of change(s):</b>
The Infrared Thermometer Daily Quality Control form now has a spot to document the date the infrared thermometer calibration expires. This is to prevent use of a thermometer with expired calibration.



- Germantown Emergency Center
- Shady Grove Medical Center
- White Oak Medical Center
- Fort Washington Medical Center

### Infrared Thermometer Daily Quality Control Form

Date Calibration Expires: \_\_\_\_\_

Month / Year: \_\_\_\_\_

Date	Test Thermometer ID	IR Thermometer ID	Test Thermometer Reading (°C)	IR Thermometer Reading (°C)	Variance (<1°C)	Acceptable? Y/N	Tech
1							
2							
3							
4							
5							
6							
7							
8							
9							
10							
11							
12							
13							
14							
15							
16							
17							
18							
19							
20							
21							
22							
23							
24							
25							
26							
27							
28							
29							
30							
31							

Reviewed by: \_\_\_\_\_

Date: \_\_\_\_\_