

 Franciscan Health System  <b>St. Clare Hospital Lakewood, WA</b> <b>St. Francis Hospital Federal Way, WA</b> <b>St. Joseph Medical Center Tacoma, WA</b> FMG Clinic Enumclaw, WA FMG Clinic Gig Harbor, WA	<h1>FORM</h1>	<b>DOCUMENT NUMBER</b> <b>R-F-UA 1038-01</b>
		Copy ID <b>J01</b>
		Unauthorized use or copying of this document is prohibited by FHS.
<b>QC/MAINTENANCE LOG FOR IRIS iQ200 MICROSCOPIC</b>		

**MONTH/YEAR:** \_\_\_\_\_

**Lot number of iQ control/Focus set** \_\_\_\_\_ **Exp. Date** \_\_\_\_\_ **Date in use** \_\_\_\_\_  
**Lot number of iQ control/Focus set** \_\_\_\_\_ **Exp. Date** \_\_\_\_\_ **Date in use** \_\_\_\_\_  
**Lot number of Lamina** \_\_\_\_\_ **Exp. Date** \_\_\_\_\_ **Date in use** \_\_\_\_\_  
**Lot number of Lamina** \_\_\_\_\_ **Exp. Date** \_\_\_\_\_ **Date in use** \_\_\_\_\_

**DAILY**

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
Check Lamina																																
New lot Lamina started																																
Run system rack																																
Perform Quality control																																
Clean surfaces																																
TECH																																

**WEEKLY**

Date				
Check Supply inventory				
TECH				

**MONTHLY**

Date	
Calibrate	
Calibrator Lot#/Exp.	
Clean Wash Bath	
Clean barcode reader/sensors	
Back up system results	
Clean in-line filter	
TECH	

**6 month:** Change peristaltic tubing(done by IRIS)

Date: \_\_\_\_\_