

## **EVOQUA WATER TESTING LOG (SJ)**

Water System: <u>ELGA MEDICA R30</u> Site: <u>4316770</u> Month/Year:							
	Out Quality Checks	RO Quality	Filter changes,	comments,			
Date	Min (>15 MΩ/cm)	(< .010 uS/cm)	and As needed r		Tech		
1							
2							
3							
4							
5							
6							
7							
8							
9							
10							
11							
12							
13							
14							
15 16							
17							
18							
19							
20							
21							
22							
23							
24							
25							
26							
27							
28							
29							
30							
31							
Monthly Check: Type 1 Water Bacterial Count (< 10 CFU/mL) Date/Tech:  Note: Notify MTC if RO Quality > .010 (ok to operate instruments up to .020).  EVOQUA - Phone: 1-800-466-7873, Account #: 04320 45716.  Comments/Reviewed by:							
J:\Lab\LAB\Document Control\Chemistry Active\Evoqua Water Testing Log (SJ)-02.doc Effective Date: 9/25/15 Page :							
Unauthorized use or copying of this document is prohibited by FHS.							

DOCUMENT APPROVAL Purpose of Document / Reason for Change:								
Changed "Note" at bottom w/ RO values. Added Evoqua phone #. New format.								
☑ No significant change to process in above revision. Per CAP, this revision does not require further Medical Director approval.								
Committee Approval Date		Medical Director Approval (Electronic Signature)	Karie Wilkinson, MD 9/25/15					