Franciscan Health System

ACL TOP - ISI VERIFICATION LOG

			dical Center Ta spital Federal		St. Clare Ho			St. Elizabeth H PSC	ospital Enum	nclaw, WA		
TERM		EXPLANATION								FORMULA		
Recombiplas	stin Lot #	Lot # of reagent currently in use on the TOP analyzer (located in analyzer settings)										
		ISI of reagent currently in use on the TOP analyzer (located in analyzer settings)									-	— ISI
ISI		ISI of reagent currently in use from package insert									Patient Pro	time
Geometric Mean		Patient mean in seconds from the TOP analyzer - called NPP: field (Normal Pool Plasma)									Seometric r	
Accession Number		Accession number of patient chosen for the INR calculation									-	
Patient Protime		Protime in seconds of patient chosen for the INR calculation										
LIS INR		INR reported in LIS on the chosen patient. Note: For INR calculation select patients with INR values >2.0 and < 5.0										
Calculated INR		INR calculated using the formula (see ACL TOP Family- ISI Weekly Verification Work Instruction)								_		
TOP INR		INR calculated by the TOP analyzer										
YEAR:												
Date	Recom	biplastin		SI Package	Geometric	Tech	Accession	Patient	LIS	Calculated	ТОР	Acceptable
MM/DD/YY	Lo	ot #	Analyzer	insert	Mean	ID	Number	Protime	INR	INR	INR	Y or N
Month:						Monthly Review By:						
Month:						Monthly Review By:						
Month:										Monthly F	Review By:	
	-											

H:\Coagulation Active\ISI Verification Log-05.docx	Effective Date: 3/15/2016	Page 1 of 1							
Unauthorized use or copying of this document is prohibited by FHS.									