



STANDARDIZED TESTING / OPERATING PROTOCOL REQUEST/ANNOUNCEMENT

New ISBT Barcodes for Special Testing of Blood Products

Description:	4 new ISBT barcodes have been defined for Blood Products labeled as Negative for Zika and CMV testing, Negative for Zika and Hemoglobin S testing, and Negative for Zika, CMV and Hemoglobin S testing. See example blood product labels.
Implementation Date:	Friday, June 23, 2017 at 11:00 am
Performing Locations:	Click on the boxes that apply: ⊠Annie Penn Hospital ⊠Moses Cone Hospital □Med Center at High Point ⊠Wesley Long Hospital ⊠Women's Hospital ⊠ Alamance Regional
Affected Locations:	Click on the boxes that apply: □ Annie Penn Hospital □ Moses Cone Hospital □ Med Center at High Point □ Wesley Long Hospital □ Women's Hospital □ Alamance Regional
Affected Departments:	Click on the boxes that apply: □ Blood Bank □ Cytology □ Flow Cytometry □ Histology □ Microbiology □ Phlebotomy □ Point of Care □ Rapid Response Lab □ Respiratory Therapy □ Specimen Processing

Specimen Type:	NA
Updated Clinical Lab Procedures:	Greensboro/Reidsville: None Alamance: None
Retired Clinical Lab Procedures:	None
Notification to Client:	Click on the boxes that apply: Section Not Applicable Memo Needed Distribution of Memo: Medical Staff Allied Health Professionals (PA, Nurse Practioners) Annie Penn (Primary Source Physicians) Dentist Emergency Department/Urgent Care Centers Family Practice Infectious Docs #ID Docs (John Campbell, Robert Comer, Jeffrey Hatcher, Cynthia Snider, Kees Van Dam) OB/GYN Pathology Pediatricians Psych Radiology Surgery #Nursing Leadership (Directors, Asst. Directors, Clinical Nurse Manager) Pharmacy - Send to DeAnne Brooks & Jim Hasspacher #IM Residents Kim Helsabeck Phlebotomy Managers and Supervisors Point of Care: Sheila, Kim & Marty
Accreditation Section:	Click on the boxes that apply: Section Not Applicable □CAP Test menu change needed □CMS Analyte form change needed □Proficiency Testing surveys changes needed or ordered

p.,	
	Click box and type needed changes/additions:
	☐Section Not Applicable
	⊠LIS changes The Sunquest consultant will make the required changes in the live
##** *#**	area at the scheduled time.
	□Reference range change/addition
Laboratory IT	☐ Technical Failure change/addition
section:	□Critical Value change/add
	☐Text comments needed
^	☐ Specimen collection instructions
	□Need to monitor TAT
	\Box CPT code for tests(s)
	On Friday, June 23 2017 at 11:00 AM, four new ISBT barcodes will be defined in SunQuest Live. The barcode labels can be scanned when labeled blood products are
	brought into blood inventory using SunQuest function Blood Product Entry (BPE).
 	The new codes should print on unit tags under attributes and should follow the unit
Technical Staff	through component processing if applicable. The new ISBT barcodes are:
Update:	1) N0127, SQ code ZIKAN: Test for Zika virus negative,
opuaic.	2) N0128, SQ code CMZIKN: CMV seronegative; Test for Zika virus negative,
	3) N0129, SQ code HSZIKN: Hemoglobin S negative; Test for Zika virus negative,
	4) N0130, SQ code CMSZKN: CMV seronegative; Hemoglobin S negative; Test for
	Zika virus negative. See barcodes in lower right corner of example blood product labels.
STOP Initiator:	Teresa Wright
	Quality Department will obtain signature:
Medical Director	flu Patimel, W 6/16/17
Signature:	plu variation, the
	Quality Department will obtain signature:
Medical Director	
Signature:	Jara (Zubinas. MD 6/20/17
	Cyarac . para de la companya del companya del companya de la compa



COMMUNITY BLOOD CENTER OF THE CAROLINAS CHARLOTTE, NC 28209 FDA Registration Number 3004054586

Properly Identity Intended Recipient
See Circular of Information for Indications,
contraindications, cautions and methods of Infusion.
This product may transmit infectious agents. Rx Only

VOLUNTEER DONOR



RED BLOOD CELLS ADENINE - SALINE (AS - 1) ADDED LEUKOCYTES REDUCED

From 500 mL CPD Whole Blood Store at 1 to 6 C









Date

RED BLOOD CELLS ADENINE – SALINE (AS – 1) ADDED LEUKOCYTES REDUCED

COMMUNITY BLOOD CENTER OF THE CAROLINAS CHARLOTTE, NC 28209

FDA Registration Number 3004054586

Properly identify intended Recipient See Circular of Information for Indications, contraindications, cautions and methods of Infusion. This product may transmit infectious agents.

Rx Only

VOLUNTEER DONOR

From 500 mL CPD Whole Blood Store at 1 to 6 C







Expiration Date



CMV NEGATIVE, ZIKA VIRUS NEGATIVE



W3985 17 062441

COMMUNITY BLOOD CENTER OF THE CAROLINAS CHARLOTTE, NC 28209 FDA Registration Number 3004054586

Properly Identify Intended Recipient See Circular of Information for Indications, contraindications, cautions and methods of Infusion. This product may transmit Infectious agents. Rx Only

VOLUNTEER DONOR



RED BLOOD CELLS ADENINE — SALINE (AS — 1) ADDED LEUKOCYTES REDUCED

From 500 mL CPD Whole Blood Store at 1 to 6 C



Rh POSITIVE

Expiration

Date



17 JUN 2017



N0130 CMV NEGATIVE, HEMOGLOBIN S NEGATIVE, ZIKA VIRUS



COMMUNITY BLOOD CENTER OF THE CAROLINAS CHARLOTTE, NC 28209

FDA Registration Number 3004054586

Properly identify intended Recipient See Circular of Information for Indications, contraindications, cautions and methods of Infusion. This product may transmit Infectious agents. Rx Only

VOLUNTEER DONOR



RED BLOOD CELLS ADENINE - SALINE (AS-1) ADDED LEUKOCYTES REDUCED

From 500 mL CPD Whole Blood Store at 1 to 6 C





Expiration Date

17 JUN 2017



HEMOGLOBIN S NEGATIVE, ZIKA VIRUS NEGATIVE