Kesler-Cribb, Sigrid

From: Kesler-Cribb, Sigrid

Sent: Tuesday, December 11, 2018 2:42 PM

To: Ballister, Kim; Barr, Kimberly; Beaman, Jeff; Bennett, Crystal (MHC-MC-Phlebotomy

> Services); Boyer, Mike; Cleveland, Thomas; Creasman, Tammy; Cruthis, Sheila; Dyson, Quanisha; Farmer, Jacee; Garcia, Amber; Hamlett, Paula; Hicks, Teri; Hogan, Donna; Holmes, Karen; Ingle, Esther (WL); James, Kelly; 'John (jdpatrick@auroradx.com)'; Kennedy, Tonya; Kish, Joshua; Kriegel, Patti; Lane, Derick; Lawson, Shelley; Lax, Amy; Lineberry, Ryan W; Manuel, Kathy; Mitchell, Shirley; Mohr, Mel; Morton, Melinda; Mullins, Tammy; O'Laughlin, Theresa; Olney, Mary S; Pegram, Jill; Roy, Sherry (Solstas); Rubinas, Tara C; Singletary, Demetria; Strickland, Marti; Turner, Wendy; Vanhoorne, Sara; Woolard, Kerri; Wright, Teresa; Younts-Meredith, Amy; Edens, Charlie; Shue, Amanda; 'Kish, Joshua'; Vestal, Mathieu; Wilkins, Vincent; Lyrla, Kathleen; Bialek, Sandy

Kesler-Cribb, Sigrid Cc:

Subject: STOP Notifications Effective December 18,2018

Attachments: 12.18.18 STOP notification Dry Ice safety ay Inspection checklist changes .pdf; 12.18.18

> STOP Notification UDS Result Comment.pdf; 12.18.18 STOP Notification ARMC moving to 24 hr DXC 660i QC.pdf; 12.18.18 STOP Notification Sunguest GUI resulting .pdf; 12.18.18 STOP Notification MC Implementing PBP2.pdf; 12.18.18 STOP Notification WL performing Lipid panel.pdf; 12.18.18 STOP Notification ARMC Xpress GAS PCR on

Infinity .pdf

STOP's Notifications Effective December 18, 2018

System Wide STOP Dry Ice safety and Safety Inspection Changes

New CAP checklist question requires a policy for the safe handling of dry ice along with documented training. The monthly and quarterly safety checklists were updated to reflect these changes along with additional changes for clarity.

System Wide STOP Drug Screen Report

Urine Drug screen report will identify the detected drugs as parent compounds, metabolites or impurities to comply with CAP updated standard.

ARMC: STOP Moving to 24 hours DXC 660i QC

ARMC main is changing from QC every 12 hours to once every 24 hours on the DxC660i.

Greensboro/Reidsville STOP Sunquest GUI resulting

Change from Flexilab OEM and MEM to GUI resulting in Sunquest.

Moses Cone Microbiology STOP PBP2 Latex Agglutination

PBP2 issued as a confirmatory test for MRSA due to an advisory notice that was sent by Biomerieux.

Wesley Long STOP WL will begin performing Lipids

Wesley Long to perform Lipid panels in house

ARMC STOP Xpress GAS (Group A Strep) PCR testing on Infinity at ARMC

Express GAS testing on the Infinity. The POCT testing rapid strep A perfume in the ED will cease and PCR testing will begin in the main laboratory.

Thanks

Sigrid© Sigrid Kesler-Cribb, MT (ASCP)

Cone Health | Laboratory Quality Manager Direct Dial: 336.832.7571 | Fax: 336.832.8270

Website: conehealth.com

		Gene	eral Information			
Title of STOP:	Drug screen report			THE REPORT OF THE PARTY OF THE	Protection was	
Requestor Name:	Kimberly Barr	nberly Barr		Submission Date:	11/30/18	
Description:	Urine drug screen report will identify the parent compounds, metabolites, or impurition: CAP updated standards		detected drugs as:	Target Date:	12/18/18	
Perforn	ning Location(s)	Affected Location(s)		Affected D		
Alamance Cancer Center High Point Cancer Center Wesley Long Cancer Center Alamance Regional Annie Penn Hospital MedCenter @ High Point MedCenter @ Mebane Moses Cone Hospital Wesley Long Hospital Women's Hospital		✓ Alamance Ca ✓ High Point C	ancer Center ancer Center Cancer Center egional Hospital High Point Mebane Hospital Hospital Hospital	Affected Department(s) Blood Bank Cytology Flow Cytometry Histology Microbiology Phlebotomy Point of Care Rapid Response Lab Respiratory Therapy Specimen Processing		
		Training / Competency				
	Training	Training Required?		YES *	✓ NO	
		Training checklist ne revised?	eds to be created or	YES *	✓ NO	
Competency Forms		Competency Log, Ob	servation form or Blind to be created or revised?	YES *	✓ NO	
	Quizzes	Quiz needs to be cre	ated or revised?	YES *	✓ NO	
* If Yes	s to the above, Explain:					
			Quality	Costion N		
	Document Control			✓ YES *	ot Applicable	
		If Yes, Have docs been submitted to MD?		YES	NO	
	* If Yes, List:	New / Revised Docs:	CHEM-0365 Beckman Ur Operating Procedure, CH Operating Procedure, CH procedure.	nicel DXC880I Integrated	el DXC600I Basic	
		Retired Documents:	NA			
		Does Current Test hav		YES *	✓ NO	
		Will New Test need a	IQCP?	YES *	✓ NO	
	* If Yes, Explain:					
N	Iotification to Client	* If yes, mark Distribu GSO/Reidsville Med ARMC Medical Staf Pharmacy - DeAnno	tion below and include modical Staff f e Brooks / Jim Hasspacher ip (Directors, Asst. Directors			

Title of STOP: Drug screen report	S spending 1100	
Memo/Notification	A STATE OF THE PARTY OF THE PAR	
Notes:		
Accreditation:	CAP Test Menu Changes Needed?	☐ YES * ✓ NO
	CMS Analyte Form Changes Needed?	YES * V NO
* 16 Van Franke	PT Survey Changes / Orders Needed?	YES * V NO
* If Yes, Explain:		

Title of STOP: Drug scre	en report					
			Billing		Section Not Applicable	
			Billing Changes Required?		YES * NO	
* If	Yes, Explain:					
		CPT/LOINC Code:				
		CDM #		00000		
		Technical Changes		CDM Modifier:		
Assigned System Adminis	trator:		mean enumges		Section Not Applicable	
Select One:	Ne Ne	ew Test / Instrumen	t	Evicting Took /	7.27	
		cate any necessary changes / additions below or mark N/A				
Reference Range:	✓ N/A	Explain Chang	es / Additions & Atta	ch patient report	example or instrument printout	
Critical Result:	✓ N/A					
Delta Check:	✓ N/A					
Technical/Clinical Range:	✓ N/A					
Autodilution:	✓ N/A					
Middleware Rules:	✓ N/A					
		A al				
Text/Comment Code:	□ IN/A	Added in "and meta	ibolites" to Urine dru	g screen text a	as applicable for reagent	
QC Code:	✓ N/A					
Specimen Requirements:	✓ N/A					
	cimen Type:	urine				
- Open	теп туре.		eterile eur			
	-		sterile cup			
Container/ Tubes: –		Acceptable:				
	-	Required Volume: 10cc				
Special Collection In	octructions.	Minimum Volume:	2cc			
Autoverification Desired:		A VEC # Total	,			
		YES * NO	N/A			
TS Tech Review been initia	ted / sched	Strains of scales are used an extended and service and	✓ YES NO	□ N/A	* Required if test is on network	
			ormation		Section Not Applicable	
Computer Systems	Involved:	Sunquest	Harvest	Order		
		PowerPath	WindowPath	Work	sheet:	
nalyzer Analyz	er Interface	YES * NO	Carried A			
nterface Required? Middleware Needed?			* If Yes, Please explain			
		YES * NO				
est Info Order Code	e Assigned?	YES * NO	Test Code(s) Assigned:			
re-Live Documentation Cor	mplete:	✓ YES	□ N/A			

Title of STOP: Drug screen report				
Medical Director Approval Quality Department will obtain signitures				
Alamance Regional Medical Director: (Signature/Date) Alamance Regional Medical Director: (Signature/Date)	MedCenter Mebane Medical Director: (Signature / Date)			
Gréensboro/Reidsville Medical Director: (Signature/Date) plu Palume, wp 13/13/18	Cancer Centers Medical Director: (Signature / Date) Not Applicable			

Approved and current. Effective starting 11/30/2018. 1719F1 (version 2.1) STOP Request/Notification Template Blank copy 2087688. Last reviewed on [n/a]. Printed on 11/30/2018 10:45 AM (EST). Page 2 of 3

STOP Notification		Name of STOP been approved by the Laboratory Medical Directors		
When:	12/18/2018			
What:		screen text will classify detected drugs as parent compounds, metabolites, or impurities.		
Who is Affected:	Location(s) Alamance Cancer Center High Point Cancer Center Wesley Long Cancer Center Alamance Regional Annie Penn Hospital MedCenter @ High Point MedCenter @ Mebane Moses Cone Hospital Wesley Long Hospital Women's Hospital	Department(s) Blood Bank Cytology Flow Cytometry Histology Microbiology Phlebotomy Point of Care Rapid Response Lab Respiratory Therapy Specimen Processing		
Why?	To comply with CAP standards updates for 2018 standards			
What you will need to do o prepare:	Nothing. New text codes will be added automatically by LIS Contact your Manager / Supervisor			
leed Help?				

Kesler-Cribb, Sigrid

From: Patrick, John <jdpatrick@auroradx.com>
Sent: Tuesday, December 11, 2018 10:07 AM

To: Kesler-Cribb, Sigrid

Subject: Re: [EXTERNAL] - please read and approve

This message was sent securely using Zix®

Sigrid,

All stops are approved. I'll sign when I'm next at MC, Thursday, 12/13 I think. Don't have schedule handy. JDP

From: Kesler-Cribb, Sigrid <Sigrid.Kesler-Cribb@conehealth.com>

Sent: Tuesday, December 11, 2018 9:50 AM

To: Patrick, John **Cc:** Kesler-Cribb, Sigrid

Subject: [EXTERNAL] - please read and approve

WARNING: This email is from an external source.

DO NOT CLICK on links or attachments unless you recognize the sender and know the content is safe. **REPORT** suspicious emails to Support@AuroraDX.com ASAP.

This message was sent securely by Cone Health.

Dr. Patrick,

I need to get STOPs out today. Can you approve these STOPs today via email and I will get you to sign the STOPs forms on Thursday when you are at Moses Cone.

Effective Date is December 18, 2018

Greensboro/Reidsville STOP Sunguest GUI resulting

Change form Flexilab OEM and MEM to GUI resulting in Sunquest

System Wide STOP Drug Screen Report

Urine Drug screen report will identify the detected drugs as parent compounds, metabolites or impurities to comply with CAP updated standard.

Moses Cone Microbiology STOP PBP2 Latex Agglutination

PBP2 issued as a confirmatory test for MRSA due to an advisory notice that was sent by Biomeriuex.

System Wide STOP Dry Ice safety and Safety Inspection Changes

Updated CAP checklist items require a policy for the safe handling of dry ice along with documented training. The monthly and quarterly safety checklist were updated to reflect these changes along with addition changes for clarity.

Wesley Long STOP WL begin performing Lipids

Wesley Long to perform Lipid panels in house

Thanks Sigrid Sigrid Kesler-Cribb, MT(ASCP)

Cone Health | Laboratory

Quality Manager

Direct Dial: 336.832.7571 | Fax: 336.832.8748

Website: conehealth.com

NOTICE: This message may contain confidential information intended only for the recipient. If you have received this communication in error, please notify the sender immediately by replying to the message and deleting it from your computer.

This message was secured by **Zix**®.

<u>CONFIDENTIALITY NOTICE</u>: This e-mail message, including any attachments, contains confidential and/or privileged information and is for the sole use of the intended recipient(s). Any unauthorized review, use, disclosure or distribution is strictly prohibited. If you are not the intended recipient and have received this e-mail in error, please contact the sender immediately by reply e-mail and destroy all copies of the original message. Thank you.

This message was secured by **Zix**®.

Greensboro/Reidsville

Example of Pre Live Report

Moses Cone Memorial Hospital 1200 N. Elm Street Greensboro, NC 27401

INTERIM REPORT

NAME: TESTING, LAB

MRN: 123456789

OPIATES

COCAINE

BARBITURATES

LOC: ED

ROOM: C29C-29 AGE: 32Y

SEX: M

{MC}

{MC}

ACCT: 111

MD: DEFAULT, PROVIDER DOB: 07/15/1986

H1866 COLL: 11/29/2018 14:15 REC: 11/29/2018 14:16 PHYS: DEFAULT, PROVIDE

DRUGS OF ABUSE SCREEN, URINE

AMPHETAMINED TETRAHYDROCANNABINOL

NONE DETECTED
NONE DETECTED [NDT] {MC} [NDT] BENZODIAZEPINES NONE DETECTED AMPHETAMINES NONE DETECTED {MC} [NDT] {MC} [NDT] {MC} NONE DETECTED [NDT]

NONE DETECTED [NDT] (NOTE)

DRUG SCREEN FOR MEDICAL PURPOSES ONLY. IF CONFIRMATION IS NEEDED FOR ANY PURPOSE, NOTIFY LAB WITHIN 5 DAYS.

LOWEST DETECTABLE LIMITS FOR URINE DRUG SCREEN

Jalus 18th no 12/7/18
Julie Patent, 2012/18/18

Drug Class Drug Class Cutoff (ng/mL)
Amphetamine and metabolites 1000 Barbiturate and metabolites 200 Benzodiazepine 200 Tricyclics and metabolites 300 Opiates and metabolites 300 Cocaine and metabolites 300 THC

 $\{MC\}$ = Performed at Moses Cone Hospital Lab, 1200 N. Elm St., Greensboro, NC 27401

TESTING, LAB

END OF REPORT PRINT DATE: 11/29/2018 PRINT TIME: 14:25

PAGE 1

Example of ARMC Pre Live Sunquest Results

Alamance Regional Medical Center INTERIM REPORT 1240 Huffman Mill Rd. Burlington, NC 27215

NAME: TESTING, LAB

AGE: 32Y SEX: M MRN: 123456789 LOC: EDA ACCT: 111 MD: DEFAULT, PROVIDER DOB: 07/15/1986

W2220 COLL: 12/12/2018 10:35 REC: 12/12/2018 10:36 PHYS: DEFAULT, PROVIDE

URINE DRUG SCREEN, QUAL			STAT
TRICYCLIC, URINE	NONE DETECTED	[NDT]	{HA}
AMPHETAMINES, URINE	NONE DETECTED	[NDT]	{AH}
MDMA, URINE	NONE DETECTED	[NDT]	{AH}
COCAINE, URINE	NONE DETECTED	[NDT]	{AH}
OPIATES, URINE	NONE DETECTED	[NDT]	{AH}
PHENCYCLIDINE, URINE	NONE DETECTED	[NDT]	{AH}
CANNABINOIDS, URINE	NONE DETECTED	[NDT]	{AH}
BARBITURATES, URINE	NONE DETECTED	[NDT]	{AH}
BENZODIAZEPINE, URINE	NONE DETECTED	[NDT]	$\{AH\}$
METHADONE, URINE	NONE DETECTED	[NDT]	{AH}
	(NOTE)		
100 Tricyclics + metabolit			
200 Amphetamines + metabol	lites, urine Cutor	ff 1000 ng/mL	

100	Tricyclics + metabolites, urine	Cutoff	1000 ng/mL
200	Amphetamines + metabolites, urine	Cutoff	1000 ng/mL
300	MDMA (Ecstasy), urine	Cutoff	500 ng/mL
400	Cocaine Metabolite, urine	Cutoff	300 ng/mL
500	Opiate + metabolites, urine	Cutoff	300 ng/mL
600	Phencyclidine (PCP), urine	Cutoff	25 ng/mL
700	Cannabinoid, urine	Cutoff	50 ng/mL
800	Barbiturates + metabolites, urine		200 ng/mL
900	Benzodiazepine, urine	Cutoff	200 ng/mL
1000	Methadone, urine	Cutoff	300 ng/mL

1200 The urine drug screen provides only a preliminary, unconfirmed 1300 analytical test result and should not be used for non-medical 1400 purposes. Clinical consideration and professional judgment should 1500 be applied to any positive drug screen result due to possible 1600 interfering substances. A more specific alternate chemical method 1700 must be used in order to obtain a confirmed analytical result. 1800 Gas chromatography / mass spectrometry (GC/MS) is the preferred 1900 confirmatory method.

{AH} = Performed at Alamance Hospital Lab, 1240 Huffman Mill Rd., Burlington, NC 27215

John B. Kul no 12/12/18-