Approved and current. Effective starting 6/28/2019. 77420.613 (version 3.1) STOP Request/Notification Template Blank copy 3474683. Last reviewed on 6/25/2019. Printed on 2/7/2020 1:26 PM (EST). Page 2 of 2

Blank copy 647 4000. Last 16 vic	:wed 011 0/25/2019. Pfilited 011 2/7/2020 1.20 PW (E	101). 1 age 2 of 2
STOP	Updated Reference Range for P2Y12	
Notification	(Platelet Inhibition) Testing	
When:	Thursday, February 13, 2020	
What?	Reference Range updated to 182-335 PRU to match current reference	
	intervals listed in the manufacturer package insert.	
Who is Affected:	Affected Location(s)	Department(s)
	Alamance Cancer Center	Blood Bank
	High Point Cancer Center	Cytology
	Wesley Long Cancer Center	Flow Cytometry
	Alamance Regional	Histology
	X Annie Penn Hospital	Microbiology
	X MedCenter @ High Point	Phlebotomy
	MedCenter @ Mebane	Point Of Care
	X Moses Cone Hospital	X Rapid Response Lab
	X Wesley Long Hospital	Respiratory Therapy
	X Women's Hospital	Specimen Processing
Why?	The initial validation reference range has been changed by the manufacturer per the latest package insert. This change is necessary to be in-line with the manufacturer.	
What you will need to do to prepare:	Review 77420.224 Reference Range List.	
	Send Post Live Documentation to IT Manager within 5 days of the effective date for test systems marked with an "X"	
	X Sunquest	WindowPath
Manager / Supervisor	X CHL	Instrument / Manual Test
Responsibility:	PowerPath	
	Manager to perform post live Patient chart review	
Need Help?	Contact your Manager / Supervisor	