


Approved and current. Effective starting 6/28/2019. 77420.613 (version 3.1) STOP Request/Notification Template  
 Blank copy 3581806. Last reviewed on 6/25/2019. Printed on 3/9/2020 10:34 AM (EDT). Page 2 of 2

 <p><b>STOP Notification</b></p>	<h2 style="text-align: center;">Fetal Fibronectin IQCP Go-Live</h2>							
<p><b>When:</b></p>	<p style="text-align: center;"><b>Tuesday, March 31, 2020</b></p>							
<p><b>What?</b></p>	<p><b>Moses Cone (MC) will implement an Individual Quality Control Plan (IQCP), effective on March 31, 2020 for Fetal Fibronectin. Wet QC will be performed with every new lot/new shipment or every 31 days, with wet QC being performed on the 15th of every month.</b></p>							
<p><b>Who is Affected:</b></p>	<p style="text-align: center;"><b>Affected Location(s)</b></p> <p>Alamance Cancer Center                  High Point Cancer Center                  Wesley Long Cancer Center                  Alamance Regional                  Annie Penn Hospital                  MedCenter @ High Point                  MedCenter @ Mebane                  X Moses Cone Hospital                  Wesley Long Hospital                  Women's Hospital</p>	<p style="text-align: center;"><b>Department(s)</b></p> <p>Blood Bank                  Cytology                  Flow Cytometry                  Histology                  Microbiology                  Phlebotomy                  Point Of Care                  X Rapid Response Lab                  Respiratory Therapy                  Specimen Processing</p>						
<p><b>Why?</b></p>	<p style="text-align: center;"><b>This is the activation date for the Fetal Fibronectin IQCP.</b></p>							
<p><b>What you will need to do to prepare:</b></p>	<p style="text-align: center;"><b>Perform wet QC with each new lot/new shipment or every 31 days, with wet QC performed on the 15th of every month.</b></p>							
<p><b>Manager / Supervisor Responsibility:</b></p>	<p style="text-align: center;">Send Post Live Documentation to IT Manager within 5 days of the effective date for test systems marked with an "X"</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%; text-align: center;">Sunquest</td> <td style="width: 50%; text-align: center;">WindowPath</td> </tr> <tr> <td style="text-align: center;">CHL</td> <td style="text-align: center;">Instrument / Manual Test</td> </tr> <tr> <td style="text-align: center;">PowerPath</td> <td style="background-color: #cccccc;"></td> </tr> </table>		Sunquest	WindowPath	CHL	Instrument / Manual Test	PowerPath	
Sunquest	WindowPath							
CHL	Instrument / Manual Test							
PowerPath								
<p><b>Need Help?</b></p>	<p style="text-align: center;">Contact your Manager / Supervisor</p>							

*This has been approved by the Laboratory Medical Directors*