Approved and current. Effective starting 7/9/2019. 77420.613 (version 3.1) STOP Request/Notification Template Blank copy 3941916. Last reviewed on [n/a]. Printed on 7/7/2020 12:30 PM (EDT). Page 2 of 2

STOP Notification	Corrections to Critical Values in SQ for CCWL and CCHP	
When:	Tuesday, July 21, 2020	
What?	Several panic values have incorrect < alerts built in SQ and do not match the CCIS	
Who is Affected:	Affected Location(s)	Department(s)
	Alamance Cancer Center	Blood Bank
	x High Point Cancer Center	Cytology
	x Wesley Long Cancer Center	Flow Cytometry
	Alamance Regional	Histology
	Annie Penn Hospital	Microbiology
	MedCenter @ High Point	Phlebotomy
	MedCenter @ Mebane	Point Of Care
	Moses Cone Hospital	x Rapid Response Lab
	Wesley Long Hospital	Respiratory Therapy
	Women's Hospital	Specimen Processing
Why?	Several panic values have incorrect < or > alerts built in SQ and do not match the CCIS. Corrections made in SQ to align with the CCIS.	
What you will need to	Creat critical value will now flag at >3.00;	
do to prepare:	Tbili critical value will now flag at >3.5;	
	ALT critical values will now flag at >275; AST critical value will now flag at >170.	
	Send Post Live Documentation to IT Manager within 5 days of the effective date for test systems marked with an "X"	
Manager / Supervisor	x Sunquest	WindowPath
Responsibility:	CHL	Instrument / Manual Test
	PowerPath	
Need Help?	Contact your Manager / Supervisor	