Form A: Wet Mount Report

Effective Date: 10/14/2019

Attach Sunquest Label:	Tech Code:
	Date/Time:
	Source:

Wet	t Mount	Finding	JS		
WBCs	None	Rare	Few	Mod	Many
Trichomonas	None	Rare	Few	Mod	Many
Yeast	None	Rare	Few	Mod	Many
Clue cells (vaginal specimens)	None	Rare	Few	Mod	Many

	Key
>10 per HPF (40x) = Many <8>	No Yeast or Trichomonas seen
5-10 per HPF (40x) = Moderate < 7 >	No Trichomonas, Clue cells, or yeast seen <n></n>
1-4 per HPF (40x) = Few <6>	No Trichomonas or Clue cells seen <v></v>
1-4 per entire slide = Rare < 5 >	Unable to r/o Trichomonas <z>- Time since collection is more than 1 hr <x></x></z>
WBCs < W >	No WBCs seen < S >
Trichomonas < T >	No Trichomonas seen <f></f>
Yeast <y></y>	No Yeast seen <h></h>
Clue cells <u></u>	No Clue cells seen <g></g>

NOTE: Always report the presence or absence of WBCs, Trichomonas, yeast, and/or Clue cells (vaginal specimens).