

## Transfusion Services Daily Temperature Log

	AM	PM	AM	PM	AM	PM	AM	PM	AM	PM	AM	PM	AM	PM
Date														
Time														

***Helmer Double Door Refrigerator (FA016406): 1-6C***

Chart														
Top														
Bottom														

***Helmer Single Door #1 (FA017203): 2-6C***

Chart														
Top														
Bottom														

***Helmer Single Door #2 (FA019668): 1-6C***

Chart														
Top														
Bottom														

***Helmer Freezer (FA029827): < -18C***

Chart														
Internal														

***Jewett Freezer (FA016405): < -18C***

Chart														
Internal														

***Platelet Rotator/Incubator (FA021411): 20-24C***

Chart														
Internal														

**NOTE:** Devices with multiple temperature readings must be in agreement within 2C. Read thermometers and chart temperatures to nearest 0.5C.

	AM	PM	AM	PM	AM	PM	AM	PM	AM	PM	AM	PM	AM	PM
Date														
Time														

**Helmer Water Bath A (FA026516): 35-37.5C**

Internal														
Digital														
H <sub>2</sub> O OK?														

**Helmer Water Bath B (FA026827): 35-37.5C**

Internal														
Digital														
H <sub>2</sub> O OK?														

**Heat Block (FA020022): 36-38C**

Position #														
Internal														

**Heat Block (FA029864): 36-38C**

Position #														
Internal														

**RT Reagent Storage Area: 1-30C**

Temp														
------	--	--	--	--	--	--	--	--	--	--	--	--	--	--

Recording tech code														
Reviewing tech code														

Comments: \_\_\_\_\_

Weekly Review: \_\_\_\_\_

Date: \_\_\_\_\_