

FORM B: HEMOCUE PLASMA/LOW HGB MAINTENANCE LOG

SERIAL # 221303011 Hemocue (B) _____

MONTH / YEAR _____

TASK	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
Daily Maintenance (after use for patient testing)																																
Clean microcuvette holder																																
Monthly Maintenance																																
Clean instrument exterior																																
As needed Maintenance																																
Clean optronic unit																																
Initials																																

CORRECTIVE ACTION	
Date/Time/Initials	Problem
	Corrective Action

Reviewed By / Date : _____

"Blank Spaces Represent Not Applicable"

Form C: Hemocue (B) Quality Control Log SN 2213303011

**QC is to be run every 8 hrs when there is a patient run, with each new lot or shipment of cuvettes, in addition to monthly. Also to be performed after a major service event, for troubleshooting.*

Microcuvettes		Lot Number	Expiration Date	Date Opened									
		<i>Check all that apply</i>											
Date	Patient Run	Monthly	New Lot or Shipment	Other	QC Level 1		QC Level 2		QC Level 3		QC OK?	Initials	Reviewed By / Date
					Lot/Exp Date:	Result	Lot/Exp Date:	Result	Lot/Exp Date:	Result			
					QC-1 Range:		QC-2 Range:		QC-3 Range:		Y / N		
					Lot/Exp Date:	Result	Lot/Exp Date:	Result	Lot/Exp Date:	Result			
					QC-1 Range:		QC-2 Range:		QC-3 Range:		Y / N		
					Lot/Exp Date:	Result	Lot/Exp Date:	Result	Lot/Exp Date:	Result			
					QC-1 Range:		QC-2 Range:		QC-3 Range:		Y / N		
					Lot/Exp Date:	Result	Lot/Exp Date:	Result	Lot/Exp Date:	Result			
					QC-1 Range:		QC-2 Range:		QC-3 Range:		Y / N		
					Lot/Exp Date:	Result	Lot/Exp Date:	Result	Lot/Exp Date:	Result			
					QC-1 Range:		QC-2 Range:		QC-3 Range:		Y / N		
					Lot/Exp Date:	Result	Lot/Exp Date:	Result	Lot/Exp Date:	Result			