

# Procedure Change/Revision (IMM 36-A)

**Procedure:** Hemoglobin A1C - Boronate Affinity HPLC Doc# IMM171

**Change/Revision:**

**Trinity has redesigned their controls and calibrators, including new catalog numbers and a shorter stability for Controls Levels 1 & Level 2.**

1. Trinity Biotech Glycated Hemoglobin Controls Level 1 & Level 2.  
Cat# 01-04-0020  
**Stability has changed from 30 days to 7 days at 2-8 °C**
2. Trinity Biotech Glycated Hemoglobin Calibrators Level 1 & Level 2  
Cat# 01-04-0022  
Stability remains at **30 days at 2-8 °C**

**Date of Revision:** 7/1/13

**Director Signature/Date:**

Handwritten signature and date: 7/5/13

**Technologist Signature/Date:** See MTS