

# Procedure Change/Revision (IMM 36-C)

**Procedure:** Quality Control Program (Doc# IMM 174)

## **Change/Revision:**

Quality Control is run daily on all assays except for batched assays where it is run once per batch. For instruments which are run for more than one shift, all QC levels are run at the beginning of first shift and a single level of QC is run at the start of second shift.

**The highlighted statement only applies if the instrument is run for a full second shift. A full second shift is defined as greater than 4 hours past the first shift. Therefore, it is not necessary to run quality control for the 12-8:30 pm shift.**

**Date of Revision:** 3/16/12

**Director Signature/Date:**

Handwritten signature of Bob Polster, dated 3/15/12.

**Technologist Signature/Date:**