

# Procedure Change/Revision (IMM 36-A)

## Procedure: Tecan Freedom EVO 100

### Change/Revision:

Changes have been made to the work table and racks have been rearranged and consolidated. Volumes to certain workgroups have been changed, as well. The chart below reflects the changes that have been made.

Work Group #	Rack #	Tube Type	Big Red Volume ( $\mu$ L)	Small Red Volume ( $\mu$ L)
Work Group 1	1	DSX, Access, ANA, etc	500	400
Work Group 2	2	Immunocap, IBT	2000	2000
Work Group 3	3	Freezer/Aliquot Tubes	500	500
Work Group 4	4	Send Outs (Mayo, Quest, ARUP)	500	500

**Date of Revision:** 1/10/12

**Director Signature/Date:**

  
3/26/2012

**Technologist Signature/Date:**