



HEMOCCULT

GVH GRAND VIEW HEALTH



HEMOCCULT SUPPLIES:

1. HEMOCCULT SENA SLIDE
2. HEMOCCULT SENA DEVELOPER
3. SPECIMEN APPLICATORS
4. GLOVES
5. GAUZE
6. TIMER
7. BIOHAZARD CONTAINER

HEMOCCULT SENSA

WHY?

HEMOCCULT SENSA IS A
DIAGNOSTIC TEST FOR RAPID
DETECTION OF FECAL OCCULT
BLOOD



HEMOCCULT RESULTS:

POSITIVE:



ANY TRACE OF BLUE ON OR AT THE EDGE OF THE SMEAR IS A POSITIVE FOR OCCULT BLOOD

Special Notes:

- There are many interfering substances such as food, drink, vitamins and medications. Therefore, this test is used for preliminary screening.
- Quality control is completed every time a test is used, if at any time the quality control check does not perform accurately, discard the slide and try using another slide and or developer.

HEMOCCULT

- ALL PATIENT TEST RESULTS ARE DOCUMENTED IN THE COMPUTER ALONG WITH THE DATE PERFORMED LOT NUMBERS AND EXPIRATION DATES. QUALITY CONTROLS ARE DOCUMENTED FOR HEMOCCULT AND GASTROCCULT TESTING.

GVH GRAND VIEW HEALTH