

Weekly Updates

May 9, 2014

Test	Description of Changes	Affected Individuals
Cerner	<p>We have had intermittent errors within Cerner and Annette @ HCMC let me know of some good tips that may prevent and/or alleviate problems.</p> <ul style="list-style-type: none">• <i>Always close open applications and then log off Cerner.</i>• <i>Then shutdown/restart/switch user on the computer.</i>	All Staff
QC on ACE Alera	<p>Just a reminder, that when you review QC results that are outside of 2 SD's, you will need to assess based on the previous result. For this, you can refer to the clipboard or check in the 3-ring binder.</p> <p>**This includes QC that is run midday when replacing reagents as well as AM QC.</p>	Bench Techs
It's that time again...	<p>Humidity has begun rising as summer approaches. That means that you will need to remove the reagent tray on Wally and wipe up the condensation with kim wipes each morning as a part of start up.</p>	Bench Techs
Flu and RSV Testing	<p>We will continue testing for FLU and RSV using the Sofia. By testing year round, we will match the process at HCMC clinics. We are working on provider education regarding year round testing.</p>	All Staff
Name Change	<p>MIALUS has changed to <u>MIALU Random</u>, which was effective Thursday May 8th.</p>	All Staff
HLA B27	<p>On May 13th there will be a pop up in EPIC alerting the providers that this can only be drawn on certain days.</p>	All Staff

If drawn on a Friday, it must be in the lab before noon. Do not draw this test on Saturday or on a Sunday of a long holiday weekend. Specimen must be received at ARUP lab within 72 hours of collection.

Procedure Additions, Updates and Changes: