

Weekly Updates

September 15, 2014

Test	Description of Changes	Affected Individuals
Update Corrected Result	I have revised this procedure. Please review the changes that are highlighted.	All Staff
Front Desk	Just a reminder: Please take YELLOW folders first. These are just verifying the patient identification and escorting to the bathroom. When necessary, please ask for help and anybody from the back can move these through the queue to keep things moving.	All Staff
Clarification: INR's	<u>INR's need to be spun and aliquotted within an hour of collection.</u> Once they are aliquotted, they are stable at RT/refrig for 24 hours.	All Staff
Millipore	A new working solution of Millipore water will be aliquotted into a new sterile cup. This has an open outdate of 7 days. Generally, this is replaced Wednesdays by opener.	Bench Techs
Dimension	Please add your initials after F/NF in the location field. This will help us track issues.	Bench Techs

Procedure Additions, Updates and Changes:

- [Corrected Results Policy](#)